




Your Reliable Dental Chair Parts Partner



Looking Forward to Working with You

 [kinsundental.com](https://www.kinsundental.com)



Foshan Kinsun Medical Instrument Co.,Ltd

- Tel: +86 18820139756
- WhatsApp: +86 18820139756
- E-Mail: Kinsundental@gmail.com
- Add: Room 418, Building D, No. 2 Lianhe Avenue, Luocun United Industrial Zone, Nanhai District, Foshan City, Guangdong Province.






About KINSUN



Foshan Kinsun Medical Instrument Co., Ltd
Your premier destination for top-quality dental chair accessories. With a comprehensive range of products including dental chair lamps, valves, suction handles, foot controls, doctor chairs, oral cameras, and dental handpiece, we are your one-stop solution for all your dental equipment needs.

At Kinsun, we pride ourselves on delivering excellence at every step. Save your time with our extensive product lineup boasting over 300 SKUs to choose from. Plus, our efficient supply chain guarantees competitive pricing without compromising on quality.



We offer quick turnaround times, with bulk orders delivered within 7 days and small orders in just 3 days. No matter the size of your order, we're committed to enhancing your efficiency and productivity.

With no minimum order quantity restrictions, complimentary samples, and OEM services available, we cater to diverse customer needs.

Whether you're a dental chair manufacturer, maintenance service provider, or second-hand equipment dealer, we've got you covered.

Partner with us today and redefine excellence in dental equipment supply.



300+ Complete Models



3-7 Days Fast Delivery



100% Quality Guaranteed



10+ Years of Experience

CONTENT

01

Air Water Syringe Series

02

Suction Handle Series

03

Dental Chair Valves

05

Dental Tubing

06

Dental Chair Plastic Parts

07

Foot Control Series

08

Dental Chair Light

09

Dental Tissue Holder

10

Doctor Chair Series

11

Oral Camera Series

12

Air Compressor Series

13

Portable Turbine Series

14

Dental Handpiece Series

15

Handheld Devices

Air Water Syringe Series

Durable, Stable and Reliable



01-01

Air-Water Syringe(Aluminum)



01-02

Air-Water Syringe(Copper)



01-03

Air-Water Syringe(SUS)



01-04

Air-Water Syringe(Green)



01-05

Air-Water Syringe(Blue)



01-06

Air-Water Syringe(Dark Blue)



01-07

Straight Syringe



01-08

Syringe Tips



01-09

Disposable Syringe Tips

Suction Handle Series

Share the Global Dentistry Market with Kinsun



02-01

Strong Suction Handle



02-02

Weak Suction Handle



02-03

Weak Suction Handle



02-04

Weak Suction Handle



02-05

Strong Suction Handle



02-06

Weak Suction Handle



02-07

Strong Suction Handle



02-08

Weak Suction Handle



02-09

Plastic Strong Suction Handle



02-10

Plastic Weak Suction Handle



02-11

Plastic Strong Suction Handle



02-12

Plastic Weak Suction Handle



02-13

Strong/Weak Suction Adaptor(Silicone)



02-14

Strong/Weak Suction Adaptor(Plastic)



02-15

Weak Suction Silicone Adaptor



02-16

Strong/Weak Suction Pipette

Dental Chair Valves

Durable, Stable and Reliable


Share the Global Dentistry Market with Kinsun



03-01 Triple Diaphragm Valve




03-02 Quadruple Diaphragm Valve



03-03 Membrane



03-04 Pressure Valve



03-05 Membrane



03-06 Pressure Valve (4 Connectors)



03-07 Pressure Valve (3 Connectors)




03-08 Regulating Valve(F)




03-09 Regulating Valve(T)




03-10 Holder Valve (open/close)




03-11 Holder Valve(close)




03-12 Holder Valve(open)



03-13 Holder Valve(close)



03-14 Holder Valve(open)



03-15 Air Switch



03-16 Scaler Water Connector



03-17 Handpiece Coupler 2 to 4



03-18 Handpiece Coupler 4 to 2



03-19 Handpiece Connector 2 Holes



03-20 Handpiece Connector 4 Holes



03-21 Electric Valve



03-22 Solenoid Valve



03-23 Main Air Switch(T)



03-24 Main Air Switch(F)



03-25 Main Air Switch(F)



04-01 Strong Suction Valve



04-02 Weak Suction Valve



04-03 Water Exchange Switch



04-04 Water Exchange Switch



04-05 Water Exchange Switch



04-06 Water-Air Switch



04-07 One Way Valve




04-08 Reducing Valve



04-09 Commutator Valve




04-10 Blowing Debris Valve




04-11 4 Holes Foot Control Valve




04-12 2 Holes Foot Control Valve



04-13 Multi-Fuction Foot Control Valve



04-14 Control Valve



04-15 Integrated Valve




04-16 Integrated Valve (18 Connectors)




04-17 Integrated Valve (left/right)



04-18 Metal Connectors



04-19 Metal Connectors



04-20 Metal Connectors



04-21 Metal Joint



04-22 Right Angle Joint



04-23 Quick Twist Connector



04-24 Quick Twist Connector























04-25 Connectors

























Dental Tubing

Durable, Stable and Reliable

Dental Chair Plastic Parts

Share the Global Dentistry Market with Kinsun

 05-01 4 Holes Handpiece Tubing	 05-02 4 Holes Handpiece Tubing (Coalescent)	 05-03 2 Holes Handpiece Tubing	 05-04 2 Holes Sprial Handpiece Tubing
 05-05 Fiber Optic Handpiece Tubing	 05-06 Syringe Tubing	 05-07 Syringe Sprial Tubing	 05-08 Strong Suction Tubing
 05-09 Weak Suction Tubing	 05-10 Suction Tubing	 05-11 4 Holes Foot Control Tubing	 05-12 4 Holes Foot Control Tubing (Coalescent)
 05-13 2 Holes Foot Control Tubing (Coalescent)	 05-14 Water / Air Tubing	 05-15 Plastic Quick Twist Joint	 05-16 Plastic Quick Twist Joint
 05-17 Plastic Quick Twist Joint	 05-18 Metal Quick Twist Joint	 05-19 Metal Quick Twist Joint	 05-20 Metal Quick Twist Joint

 06-01 Water Heater	 06-02 Water Heater	 06-03 Square Pressure Meter	 06-04 Round Pressure Meter (Big/Small)
 06-05 Fiber Optic Box	 06-06 Fiber Optic Box (Taiwan)	 06-07 Air Reducer (Round)	 06-08 Air Reducer
 06-09 Air Reducer (Square)	 06-10 Transformer (80w/200w)	 06-11 Water Filter (Copper)	 06-12 Water Filter (Platic)
 06-13 Strong/Weak Suction Filter	 06-14 Strong and Weak Suction Filter	 06-15 Strong Suction Sliencer	 06-16 Flush/Supply Water Tube
 06-17 Spittoon Filter	 06-18 Spittoon Filter (4pcs/set)	 06-19 Spittoon Filter (5pcs/set)	 06-20 Plastic Adaptor
 06-21 Plastic Adaptor	 06-22 Storage Bottle	 06-23 Bottle Cap	 06-24 Air Bottle

Foot Control Series

Durable, Stable and Reliable



07-01
2 Holes Foot Control



07-02
4 Holes Foot Control



07-03
Multi-Fuction Foot Control



07-04
Multi-Fuction Foot Control



07-05
Multi-Fuction Foot Control



07-06
Square Foot Control



07-07
Multi-Fuction Foot Control



07-08
Multi-Fuction Foot Control



07-09
Foot Control Switch



07-10
Dental Chair Motor



07-11
Dental Chair Motor



07-12
Dental Chair Control System

Dental Chair Light

Share the Global Dentistry Market with Kinsun



08-01
Dental Halogen Light



08-02
LED Reflector Light



08-03
Dental LED Light



08-04
Dental LED Light



08-05
Dental LED Light



08-06
Dental LED Light



08-07
Dental LED Light



08-08
Implant LED Light



08-09
Light Arm



08-12
Film Viewer



08-13
Ponorama Film Viewer



08-10
Light Arm



08-14
Wireless Bulb



08-15
Bulb



08-11
Light Arm



08-16
Glass Reflector
(Square/Round)



08-17
Lamp Switch

Dental Tissue Holder

Durable, Stable and Reliable



09-01
3 in 1 Tissue Holder(Deluxe)



09-02
3 in 1 Tissue Holder(Blue)



09-03
3 in 1 Tissue Holder(White)



09-04
Tissue Box(Blue/White)



09-05
Blue Cup Cylinder




09-06
White Cup Cylinder




09-07
Square Tray



09-08
Arc Tray



09-09
Handpiece Holder(Big)



09-10
Handpiece Holder(Small)




09-11
Single Holder(Big)



09-12
Single Holder(Small)



09-13
Glass Spittoon



09-14
Glass Spittoon



09-15
Ceramic Spittoon



09-16
Glass Spittoon

Doctor Chair Series

Share the Global Dentistry Market with Kinsun



10-01
Doctor Chair



10-02
Doctor Chair



10-03
Doctor Chair



10-04
Doctor Chair



10-05
Doctor Chair




10-06
Doctor Chair



10-07
Doctor Chair




10-08
Doctor Chair




10-09
Doctor Chair Cushion



10-10
Gas Spring



10-11
Chair Chassis



10-12
Chair Chassis



10-13
Plastic Chair Foot



10-14
Metal Chair Foot



10-15
Chair Wheels



10-16
Chair Wheels

Oral Camera Series

Durable, Stable and Reliable



11-01
Oral Camera

HD Oral Camera

Screen: 1080P
 Pixels: 5 Mega Pixels
 Focus Distance: 3-50mm Auto Focus
 Image Sensor: 1/4 SONY HAD CCD
 Camera: 6 LED Handpiece
 NTSC/PALSigna System
 Wifi Fuction Available
 8GB USB Storage
 Dynamic Anti Shake Shooting



11-02
Oral Camera

Computer Oral Camera System

Screen: 24 Inch, 1920*1080, Touch Screen
 CPU: Intel Core i3/2.40 Ghz
 Storage: SSD 256GB
 Ram: 8GB
 System: Windows 10, 64bt
 Camera: 1500 Mega Pixels
 Adapter: USB*4, VGA*1, HDMI*1
 Wifi Fuction Available
 Can be connected to X-ray sensors, CT, etc.



11-03
Oral Camera Holder



11-04
Oral Camera Holder



11-05
Oral Camera Holder

Air Compressor Series

Share the Global Dentistry Market with Kinsun



12-01
Air Compressor



12-02
Air Compressor



12-03
Air Compressor



12-04
Air Compressor

Model	Air Tank	Power	Air Flow	For Dental Chair	Weight	Noise	Packing Size	Voltage
12-01	32L	600W	70L/min	1 Unit	25KG	55-62dB	45*45*65cm	110V/220V(50/60Hz)
12-02	35L	850W	100L/min	1-2 Unit	26KG	55-64dB	45*45*67cm	110V/220V(50/60Hz)
12-03	60L	1100W	140L/min	2-3 Unit	45KG	60-65dB	45*45*85cm	110V/220V(50/60Hz)
12-04	60L	1700W	200L/min	3-4 Unit	47KG	60-68dB	45*45*85cm	110V/220V(50/60Hz)



12-04
Mobile Suction Unit

Voltage: 220v/110v; 50Hz/60Hz
 Power: 550w
 Fan Speed: 2800r/min
 Noise: 53dB
 Vacuum: -11Kpa
 Air Flow: 100L/min
 Weight: 35KG
 Size: 37*33*83cm
 Package:45*43*96cm



12-05
Mobile Suction Unit

Voltage: 220v/110v; 50Hz/60Hz
 Power: 250w
 Noise: ~13KPa
 Control: Touch Screen
 Gross: 50 Kg
 Size: 52*53*90cm
 2 Tips+1 Adapter
 Removable Filter
 Auto Cleaning

Portable Tubine Series

Durable, Stable and Reliable




13-01
Portable Chiar

- Configuration:
- Folding Chair * 1
- Plastic Spittoon * 1
- LED light * 1
- Metal tray * 1
- Size:97*64.5*32.5



13-02
Portable Chiar

- Configuration:
- Folding Chair * 1
- Plastic Spittoon * 1
- LED light * 1
- Metal tray * 1
- Turbine Machine * 1



Color for Portable Chair

Medical PU Leather
Healthy and Comfortable
Folding chairs come in multiple colors
to meet your personalized needs
Note: Special colors may have a
minimum order quantity requirement



13-03
Portable Turbine Machine



13-04
Portable Turbine Machine



13-05
Mobile Turbine Machine



13-06
Turbine Machine



13-07
Turbine Machine



13-08
Turbine Machine(Metal)



13-09
Turbine Machine(Platic)

Dental Handpiece Series

Share the Global Dentistry Market with Kinsun



14-01

High Speed Handpiece

- Single Water Spray
- Ceramic Ball Bearing
- Standard Head
- Push Button
- 2/4-Hole Optional
- 135°C Disinfection
- Without Led Type



14-02

High Speed Handpiece

- Single Water Spray
- Ceramic Ball Bearing
- Standard Body
- Push Button
- 2/4-Hole Optional
- 135°C Disinfection
- Without Led Type



14-03

High Speed Handpiece

- Triple Water Spray
- Ceramic Ball Bearing
- Stainless Body
- Push Button
- 2/4-Hole Optional
- 135°C Disinfection
- Led Type



14-04

High Speed Handpiece

- 5 Water Spray
- Ceramic Ball Bearing
- Shadowless 5 Led Lamp
- Push Button
- 2/4-Hole Optional
- 135°C Disinfection
- Non-Slip Handle



14-05

Low Speed Handpiece

- Outer Water Spray
- Ceramic Ball Bearing
- Super Short Motor
- Push Button
- 1:1 Ratio
- 2/4-Hole Optional
- 135° C Disinfection



14-06

Low Speed Handpiece

- Outer Water Spray
- Ceramic Ball Bearing
- Long Air Motor
- Screw Type
- 1:1 Ratio
- 2/4-Hole Optional
- 135° C Disinfection



14-07

Implant Handpiece (Optic)

- 20:1 Ratio
- Max Speed: 2000rpm
- Torque: 80 Ncm
- Titanium Body
- Pushn Button Chuck
- Double Sealing System
- Cellular Glass Optics



14-08

Implant Handpiece

- 20:1 Ratio
- Max Speed: 2000rpm
- Torque: 80 Ncm
- Stainless Steel Body
- Pushn Button Chuck
- Double Sealing System
- No Lighting System

Handheld Devices

Durable, Stable and Reliable



15-01

LED Curing Light

Wave Length : 385-515nm
 Mode 1: 1200mW/cm²
 Mode 2: 2300mW/cm²
 Color: 2 Colors



15-02

LED Curing Light

Wave Length: 420-480nm
 Mode 1: 1200mW/cm²
 Mode 2: 1500mW/cm²
 Color: 5 Colors



15-03

Gutta Obturation Heated Pen

Power: 10W
 Charging time: 4 Hours
 Battery: 1.2V 2400mAh
 Operating Voltage: DC.3.6V



15-04

Gutta Cutter

Heat Time: 2s
 Max Temp: 260° C
 Battery: 2400mAh
 Tips: 4pcs



15-05

Sandblast Machine

Stable, Reliable, Save Time
 Easy to Operate
 Easy to Maintain
 Quickly Connect



15-06

Ultrasonic Endo Activate

Working Frequency: 45KHz
 Battery: 1500mAh(DC3.7V)
 Atmospheric Pressure: 70kPa ~106kPa
 Contra Angle Design for Increased Accessibility



15-07

Light Guide Rod



15-08

Gutta Cutter Tips



15-09

Scaler Tips

FAQS

Share the Global Dentistry Market with Kinsun

Q: Whats the minimum order quantity of your products?

A: Our dental products can be purchased in small quantities, but you should also consider the cost of shipping. Generally speaking, buying in large quantities can save a lot.

Q: Can you provide free samples?

A: Yes, we know that samples are needed for testing before confirming large orders, so some of our products can provide free samples, but you need to bear the shipping cost yourself.

Q: Can you do dental products OEM/ODM?

A: Yes, we can provide you OEM/ODM of products, We can do logo, packaging, stickers, or product development.

Q: How long is the delivery time?

A: The delivery time for most products is 3-10 days, if you have a large order and a large variety of products, it may take up to 15-20 days.

Q: How long is your products warranty?

A: The warranty period of each product is different, please contact us for specific details.

Q: What payment methods do you accept?

- A: 1. US dollar company account
 2. US dollar personal account
 3. Western Union
 4. RMB company account
 5. RMB Private Account
 6. Alibaba pay